




Specimen Packaging Instructions



- (1) Please complete:
- Requisition
 - Insurance Information
 - Specimen bottle labels

Label and place specimen container in provided biohazard bag.

(Optional) Remove yellow copy of requisition for your records.



- (2) Please construct included white shipping box:
- Orient box as shown
 - Fold up bases at the bottom of the box, then fold in 2 tabs located at each end of the box
 - Fold up the 2 remaining long flaps, and tuck them into the 2 slits located at the bottom of the box
 - Fold in the 2 curved flaps to close the lid of the box



- (3) Fold requisition form and insert into the pocket of the biohazard bag.

Place biohazard bag(s) into the box, close the box, and place box(es) into FedEx Clinical Pak.

Mailing:

- Seal FedEx Clinical Pak**
- Attach FedEx Airbill to package
- Call FedEx for pickup or Drop off at FedEx location